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#### ●Motor Change

Standard motor can be changed to customer specified motor by request.

\*Depends on the exchanged motor's specification, that of specification way not be achieved the catalog specification.

\*If customer specified motor is not stepping motor, that of inspection will be done by using standard motor only.

#### ●Change Price

①Same size motor change price:In case of same size motor's mounting and shaft diameter as standard motor's.

②Different size motor change price:In case of different size motor's mounting or shaft diameter from standard motor's.

\*Change price does not include the motor price. Thus additional motor price would be charged without supplied motor.

\*Different motor price estimated with no motor's disc (Home Position Sensor) type, in case of product with disc (Home position sensor) type. If it is necessary, additional cost would be charged.

(In this case, Proximity home position sensor equipped with motorized stage will be performed as home position sensor.)

\*Depending on the motor's type, motor price will be changed. Please inform us the motor model.

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
### Vacuum category and pressure range

	区分	圧力範囲 (Pa)	圧力範囲 (Torr)
低真空	Low Vacuum	100kPa~100Pa	760~1 Torr
中真空	Medium Vacuum	100Pa~0.1Pa	1~10 <sup>-3</sup> Torr
高真空	High Vacuum	0.1Pa~10 <sup>-3</sup> Pa	10 <sup>-3</sup> Torr~10 <sup>-7</sup> Torr
超高真空	Ultra-High Vacuum	10 <sup>-5</sup> Pa~10 <sup>-8</sup> Pa	10 <sup>-7</sup> Torr~10 <sup>-10</sup> Torr
極高真空	Extreme-High Vacuum	≤10 <sup>-8</sup> Pa	<10 <sup>-10</sup> Torr

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When you mount the stage onto other stage or a bread board, you need to connect the stage to a motor controller or a motor driver, then move the stage's table until appear the mounting hole for the screw. If that of mounting hole position does not match your favor, would you check the translation spacer in various size (see page O-006).

\*Depend on the stage, the mounting hole does not appear only using motor driving. In that case, move the stage with motor-shaft rotating by manual, but there is possibility of heavy damage like sensor destruction caused by limit over moving. When the mounting hole is able to use by manual moving, please note to operate carefully by manual.

Stage lubricants can be changed for clean-room and vacuum applications.

\*Except lubricant, using parts and surface treatment and manufacturing environment is same as standard condition.

\*Clean room environment or vacuum environment cannot be guaranteed by option lubricant change.

### ●Clean room lubricant: LG2 from NSK LTD.

(Clean room lubricant aids in the reduction of particulate output when stage's screw and bearing rotating.)

### ●Vacuum lubricant: S-3230 from MORESCO corporation.

(Vacuum oil is low in particle emissions, and Low saturated vapor pressure. The oil is available, when you don't care the contamination or the case of Helium displacement by rotary pump suction. More clean and high vacuum stages are lined up. Refer to page M-002 for more information.)

\*Viscosity is different from standard lubricant, then feeling can be different.

\*Customer's specified lubricant can be applied (lubricant cost will be added separately).

### ●Affected parts

1. Guide: ball race, sliding surface of ball bearing, linear guide, crossed roller bearing.
2. Feeding mechanism: ground screw, ball screw, worm and worm wheel, wedge, rack and pinion
3. Bearing

\*Inside of motor grease can not be changed.

\*About the linear guide or ball screw, grease change will be done after simplified cleaning in kohzu.

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### conversion chart

	(Torr)	(atm)	(kgf/cm <sup>2</sup> )	(psi)	(Pa)
1Torr=1mmHg	1	1.31579×10 <sup>-3</sup>	1.35951×10 <sup>-3</sup>	1.9337×10 <sup>-2</sup>	1.33322×10 <sup>2</sup>
1atm	760	1	1.03323	1.4695×10	1.01325×10 <sup>5</sup>
1kgf/cm <sup>2</sup>	7.3556×10 <sup>2</sup>	9.6784×10 <sup>-1</sup>	1	1.4223×10	9.80665×10 <sup>4</sup>
1psi	5.1715×10	6.8046×10 <sup>-2</sup>	7.0307×10 <sup>-2</sup>	1	6.8948×10 <sup>3</sup>
1Pa (=1N/m <sup>2</sup> )	7.5006×10 <sup>-3</sup>	9.8692×10 <sup>-6</sup>	1.01972×10 <sup>-5</sup>	1.4503×10 <sup>-4</sup>	1